

Product Datasheet

Mouse UCHL1 ELISA Kit (orb442815)

Description	Mouse Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9R0P9
Sensitivity	25pg/mL
Entrez	22223
Expiration Date	Please enquire.