

## Product Datasheet

### Rat UCHL1 ELISA Kit (orb442816)

<b>Description</b>	Rat Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) Sandwich ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	78.125-5000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q00981</a>
<b>Sensitivity</b>	25pg/mL
<b>Entrez</b>	<a href="#">29545</a>
<b>Expiration Date</b>	Please enquire.