

## Product Datasheet

### General Vitamin A ELISA Kit (orb442864)

<b>Description</b>	General Vitamin A (VA) competitive inhibition ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Sensitivity</b>	0.03ng/mL
<b>Expiration Date</b>	Please enquire.