


Product Datasheet

Human Visfatin ELISA Kit (orb442894)

Description	Human Visfatin (VF) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P43490
Sensitivity	6.3pg/mL
Entrez	10135
Expiration Date	Please enquire.



Standard
Curve of
Human
Visfatin
ELISA K...