

Product Datasheet

Human VIP ELISA Kit (orb442905)

| | |
|--------------------------|---|
| Description | Human Vasoactive Intestinal Peptide (VIP) competitive inhibition ELISA Kit |
| Reactivity | Human |
| Range | 7.812-500pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | P01282 |
| Sensitivity | 2.34pg/mL |
| Entrez | 7432 |
| Expiration Date | Please enquire. |