

## Product Datasheet

### Mouse VIP ELISA Kit (orb442906)

<b>Description</b>	Mouse Vasoactive Intestinal Peptide (VIP) competitive inhibition ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	7.812-500pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>P32648</b>
<b>Sensitivity</b>	2.26pg/mL
<b>Entrez</b>	<b>22353</b>
<b>Expiration Date</b>	Please enquire.