

Product Datasheet

Rat VIP ELISA Kit (orb442908)

Description	Rat Vasoactive Intestinal Peptide (VIP) competitive inhibition ELISA Kit
Reactivity	Rat
Range	7.812-500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01283
Sensitivity	2.64pg/mL
Entrez	117064
Expiration Date	Please enquire.