

Product Datasheet

Human VSIG2 ELISA Kit (orb442921)

Description	Human V-Set And Immunoglobulin Domain Containing Protein 2 (VSIG2) competitive inhibition ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q96IQ7
Sensitivity	0.152ng/mL
Entrez	23584
Expiration Date	Please enquire.