

Product Datasheet

Mouse vWF ELISA Kit (orb442933)

Description	Mouse Von Willebrand Factor (vWF) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	0.25-16ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q8CIZ8
Sensitivity	0.11ng/mL
Entrez	22371
Expiration Date	Please enquire.