

## Product Datasheet

### Porcine vWF ELISA Kit (orb442934)

<b>Description</b>	Porcine Von Willebrand Factor (vWF) competitive inhibition ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	15.625-1000ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	plasma and other biological fluids.
<b>Sample Types</b>	plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q28833</a>
<b>Sensitivity</b>	4.51ng/mL
<b>Entrez</b>	<a href="#">399543</a>
<b>Expiration Date</b>	Please enquire.