

Product Datasheet

Human WNT5B ELISA Kit (orb442964)

Description	Human Wingless Type MMTV Integration Site Family, Member 5B (WNT5B) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9H1J7
Sensitivity	0.53ng/mL
Entrez	81029
Expiration Date	Please enquire.