

Product Datasheet

Human ZFAND6 ELISA Kit (orb442994)

Description	Human Zinc Finger, AN1-Type Domain Protein 6 (ZFAND6) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6FIF0
Sensitivity	0.056ng/mL
Entrez	54469
Expiration Date	Please enquire.