

Product Datasheet

Human ZG16B ELISA Kit (orb442997)

Description	Human Zymogen Granule Protein 16 Homolog B (ZG16B) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96DA0
Sensitivity	0.131ng/mL
Entrez	124220
Expiration Date	Please enquire.