

Product Datasheet

Chicken Glucagon ELISA Kit (orb443003)

Description	Chicken Glucagon (GC) competitive inhibition ELISA Kit
Reactivity	Gallus
Range	23.437-1500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P68259
Sensitivity	7.25pg/mL
Entrez	396196
Expiration Date	Please enquire.