



Product Datasheet

HMGN2 Antibody (orb443184)


Description	HMGN2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human HMGN2 (KKGKADAGKEGNNPAENGDAKTDQAQKAEGAGDAK).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	P05204
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of Raji cells us...



WB analysis of HMGN2 using anti-HMGN2 an...



IF analysis of HMGN2 using anti-HMGN2 an...