



Product Datasheet

LATS1 Antibody (orb443195)


Description	LATS1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC, IHC-Fr, WB
Immunogen	E. coli-derived human LATS1 recombinant protein (Position: Q637-A698).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	150 kDa
Uniprot ID	O95835
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of SiHa cells using anti-LATS1 antibody.



Western blot analysis of LATS1 using anti-LATS1 antibody.



IHC analysis of LATS1 using anti-LATS1 antibody.