



Product Datasheet

Rad51 Antibody (orb443214)


Description	Rad51 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human Rad51 (KKLEEAGFHTVEAVAYAPKKELINIKGISEAKADK).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q06609
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of SiHa cells us...



WB analysis of Rad51 using anti-Rad51 an...



WB analysis of Rad51 using anti-Rad51 an...