


Product Datasheet

Human AGT/Serpin A8/Angiotensinogen ELISA Kit (orb443235)

| | |
|----------------------------|---|
| Description | Human AGT/Serpin A8/Angiotensinogen ELISA Kit |
| Reactivity | Human |
| Range | 312 pg/ml - 20,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: D34-A485 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:500, concentration in serum and plasma is 1-2 ug/ml |
| Sample Types | cell culture supernatants, serum and plasma (EDTA). |
| Uniprot ID | P01019 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Human AGT PicoKine ELISA Kit standard CU...