


Product Datasheet

Human NT-proANP/NPPA ELISA Kit (orb443239)

Description	Human NT-proANP/NPPA ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: N26-R123
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is 0.77-45.96 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P01160
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human NT-proANP/NPPA PicoKine ELISA Kit ...