


Product Datasheet

Mouse TRKC ELISA Kit (orb443251)

Description	Mouse TRKC ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C32-T429
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is around 300pg/ml
Sample Types	cell culture supernatants and cell lysates.
Uniprot ID	Q6VNS1
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse TRKC
PicoKine
ELISA Kit
standard c...