


Product Datasheet

Mouse CRLF2/TSLP R ELISA Kit (orb443256)

Description	Mouse CRLF2/TSLP R ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A20-L233
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 400 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q8CII9
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CRLF2/TSLP
R PicoKine
ELISA Kit
st...