


Product Datasheet

Human CDCP1 ELISA Kit (orb443258)

Description	Human CDCP1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F30-T667
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
Uniprot ID	Q9H5V8
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human
CDCP1
PicoKine
ELISA Kit
standard ...