

Product Datasheet

PRKD1 (phospho-Ser738) antibody (orb4438)

| | |
|----------------------------|--|
| Description | PRKD1 (phospho-Ser738) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human PRKD1 around the phosphorylation site of Ser738 (EK(p-S)FR) |
| Target | phospho-PRKD1 (Ser738) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 102 kDa |
| Uniprot ID | Q15139 |
| Dilution Range | IHC-P=1:100-500 |
| Expiration Date | 12 months from date of receipt. |