



Product Datasheet

HLA-G antibody (orb44431)


Description	Mouse Monoclonal to HLA-G.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Bacterially expressed recombinant HLA-G heavy chain (denatured).
Target	HLA-G
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution 1-2 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P17693
Entrez	3135
Dilution Range	Western blotting: Recommended dilution 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Western blotting analysis of HLA-G in LC...



Immunohistochemistry staining of kidney ...



Western blotting analysis of human HLA-G...