

Product Datasheet

HIV-1 gp24 antibody (orb44463)

Description	Mouse Monoclonal to HIV-1 gp24.
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	HIV-1 p24 C-terminal peptide.
Target	Capsid protein p24 (HIV-1)
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Expiration Date	12 months from date of receipt.