


## Product Datasheet

### STAT1 (phospho-Ser727) antibody (orb44474)

<b>Description</b>	Mouse Monoclonal to STAT1 (phospho-Ser727).
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP, WB
<b>Immunogen</b>	STAT1 peptide sequence 721-733 (Ser727 phosphorylated).
<b>Target</b>	STAT1 (pS727)
<b>Preservatives</b>	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 2-8°C. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified by protein-A affinity chromatography.
<b>Uniprot ID</b>	<a href="#">P42224</a>
<b>Entrez</b>	<a href="#">6772</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blotting analysis (reducing cond...)