


Product Datasheet

GST antibody (orb44500)

Description	Mouse Monoclonal to GST.
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	GST fusion protein
Target	GST tag
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 1 µg/ml; positive control: GST fusion protein, reducing conditions.
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Western blotting: Recommended dilution: 1 µg/ml; positive control: GST fusion protein, reducing conditions.
Expiration Date	12 months from date of receipt.



Western blotting analysis of isolation o...