


Product Datasheet

HSV1 + HSV2 antibody (orb44502)

| | |
|----------------------------|---|
| Description | Mouse monoclonal to HSV1+HSV2. |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, ICC |
| Immunogen | Extract of HSV-1 infected VERO (green monkey kidney) cells. |
| Target | Glycoprotein B (HSV-1/HSV-2) |
| Preservatives | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide |
| Concentration | 1 mg/ml |
| Storage | Store at 2-8°C. Do not freeze. |
| Note | For research use only |
| Isotype | Mouse IgG2b |
| Clonality | Monoclonal |
| Purity | Purified by protein-A affinity chromatography. |
| Uniprot ID | P10211 |
| Entrez | 24271469 |
| Expiration Date | 12 months from date of receipt. |



Immunocytochemistry detection of HSV by ...