

## Product Datasheet

### HSV1 antibody (orb44505)

<b>Description</b>	Mouse Monoclonal to HSV1.
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, WB
<b>Immunogen</b>	Extract of HSV-1 infected VERO (green monkey kidney) cells.
<b>Target</b>	Glycoprotein C (HSV-1)
<b>Preservatives</b>	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 2-8°C. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
<b>Uniprot ID</b>	<a href="#">Q8UYE2</a>
<b>Entrez</b>	<a href="#">2703410</a>
<b>Expiration Date</b>	12 months from date of receipt.