




Product Datasheet

gamma Tubulin antibody (orb44552)

Description	Mouse Monoclonal to gamma Tubulin.
Reactivity	Animal, Plant, Protozoa
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	human gamma-tubulin peptide EYHAATRPDYISWGTO, amino acids 434-449
Target	gamma-Tubulin
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 1-2 µg/ml, reducing conditions. Immunocytochemistry: Methanol/acetone fixation required.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P23258
Entrez	7283
Dilution Range	Western blotting: Recommended dilution: 1-2 µg/ml, reducing conditions. Immunocytochemistry: Methanol/acetone fixation required.
Expiration Date	12 months from date of receipt.



Western blotting analysis of differentia...



Western blotting analysis of human gamma...