


## Product Datasheet

### CD16 antibody (Biotin) (orb44654)

<b>Description</b>	Mouse Monoclonal to CD16.
<b>Reactivity</b>	Human, Primate
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	FC, IHC-Fr, IP
<b>Immunogen</b>	Human neutrophils
<b>Target</b>	CD16
<b>Preservatives</b>	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 2-8°C. Do not freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	Flow cytometry: Recommended dilution: 5 µg/ml.
<b>Isotype</b>	Mouse IgG1 kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
<b>Dilution Range</b>	Flow cytometry: Recommended dilution: 5 µg/ml.
<b>Expiration Date</b>	12 months from date of receipt.



Flow cytometry analysis (surface stainin...