



Product Datasheet

CD19 antibody (orb44660)

Description	Mouse Monoclonal to CD19.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	Human CCL (chronic lymphocytic leukemia) cells
Target	CD19
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P15391
Entrez	930
Dilution Range	Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD19 positive lympho...



Flow cytometry surface staining pattern ...