

Product Datasheet

ANGPTL4 antibody (Biotin) (orb446828)

Description	ANGPTL4 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ANGPTL4/ARP4
Target	Angiopoietin-like 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BY76
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.