



Product Datasheet

CD229 antibody (orb44727)

Description	Mouse Monoclonal to CD229.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	CD299-transfected 300.19 pre-B cell line
Target	CD229
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 6 µg/ml; positive control: peripheral blood lymphocytes.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9HBG7
Entrez	4063
Dilution Range	Flow cytometry: Recommended dilution: 6 µg/ml; positive control: peripheral blood lymphocytes.
Expiration Date	12 months from date of receipt.



Separation of human CD229 positive lymph...



Flow cytometry surface staining pattern ...