

Product Datasheet

EGFP Tag antibody (Biotin) (orb447439)

Description	EGFP Tag antibody (Biotin)
Species/Host	Rabbit
Reactivity	All
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant EGFP protin
Target	EGFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at -80 °C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42212
Dilution Range	WB=1:1000-10000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.