

Product Datasheet

Glu-Glu Tag antibody (Biotin) (orb447631)

Description	Glu-Glu Tag antibody (Biotin)
Species/Host	Rabbit
Reactivity	All
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated Glu-Glu Tag
Target	Glu-Glu Tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:1000-5000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.