

Product Datasheet

HA Tag(12CA5) antibody (Biotin) (orb447732)

Description	HA Tag(12CA5) antibody (Biotin)
Species/Host	Mouse
Reactivity	All
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated HA tag
Target	HA tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB=1:5000-50000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.