



Product Datasheet

Ly6G antibody (FITC) (orb44841)

Description	Rat Monoclonal to Ly6G.
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Murine granulocytes
Target	Ly6G
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Rat IgG2b
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	A0A087WQF5
Entrez	546644
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of murine Ly6G positive cells...



Flow cytometry surface staining pattern ...