

Product Datasheet

PSBO antibody (HRP) (orb44844)

Description	PSBO antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Spinacia oleracea Oxygen-evolving enhancer protein 1, chloroplastic protein (63-232AA)
Target	PSBO
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P12359
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.