

## Product Datasheet

### PSBO antibody (FITC) (orb44845)

|                        |  |
|------------------------|--|
| <b>Description</b>     | PSBO antibody (FITC)   |
| <b>Species/Host</b>    | Rabbit   |
| <b>Reactivity</b>      | Plant  |
| <b>Conjugation</b>     | FITC   |
| <b>Immunogen</b>       | Recombinant Spinacia oleracea Oxygen-evolving enhancer protein 1, chloroplastic protein (63-232AA) |
| <b>Target</b>          | PSBO   |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol  |
| <b>Form/Appearance</b> | Liquid   |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                      |
| <b>Note</b>            | For research use only  |
| <b>Isotype</b>         | IgG  |
| <b>Clonality</b>       | Polyclonal   |
| <b>Uniprot ID</b>      | <a href="#">P12359</a>   |
| <b>Dilution Range</b>  | ELISA  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |