

Product Datasheet

VACWR062 antibody (Biotin) (orb44865)

Description	VACWR062 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Vaccinia virus Protein E6 protein (180-429AA)
Target	VACWR062
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21607
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.