

## Product Datasheet

### YES1 (phospho-Tyr537) antibody (orb4498)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | YES1 (phospho-Tyr537) antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, IF, IHC-Fr, IHC-P, WB  |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human YES1 around the phosphorylation site of Tyr537 (PQ(p-Y)QP) |
| <b>Target</b>              | phospho-YES1 (Tyr537)   |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.  |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 60 kDa  |
| <b>Uniprot ID</b>          | <a href="#">P07947</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |