

Product Datasheet

LENG4 antibody (Biotin) (orb451165)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody against LENG4 conjugated to Biotin |
| Species/Host | Rabbit |
| Conjugation | Biotin |
| Tested Applications | ELISA, ICC, IHC-Fr, IHC-P |
| Immunogen | KLH conjugated synthetic peptide derived from human LENG4 |
| Target | MBOAT7 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Lyophilized or Liquid |
| Concentration | - 100 µg (in 100 µl): 1 mg/ml |
| Storage | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 53kDa |
| Uniprot ID | Q96N66 |
| Dilution Range | IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |