

Product Datasheet

Tubulin Alpha antibody (Biotin) (orb452727)

| | |
|----------------------------|--|
| Description | Tubulin Alpha antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Biotin |
| Tested Applications | ELISA, ICC, IHC-Fr, IHC-P, WB |
| Immunogen | KLH conjugated synthetic peptide derived from human TUBA1A + TUBA1B |
| Target | Tubulin alpha (Loading Control) |
| Preservatives | Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q71U36 |
| Dilution Range | WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |