

Product Datasheet

Rabbit IgG(Fc Specific,Heavy Chain) antibody (Biotin) (orb453045)

Description	Rabbit IgG(Fc Specific,Heavy Chain) antibody (Biotin)
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	rabbit IgG
Target	Rabbit IgG(Fc specific,Heavy Chain)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB=1:500-1000,IHC-P=1:100-500,IHC-F=1:100-500
Expiration Date	12 months from date of receipt.