

Product Datasheet

Alpha Tubulin antibody (Biotin) (orb453347)

Description	Alpha Tubulin antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant human alpha Tubulin Protein
Target	alpha Tubulin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P68366
Dilution Range	WB=1:50000-200000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.