

Product Datasheet

Alpha Tubulin (Acetyl Lys40) antibody (Biotin) (orb453348)

Description	Alpha Tubulin (Acetyl Lys40) antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Mouse, Zebrafish
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised acetylpeptide derived from human alpha Tubulin around the acetylation site of Acetyl Lys40
Target	alpha Tubulin (Acetyl Lys40)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P68366
Dilution Range	WB=1:500-1000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.