

Product Datasheet

Maltose Binding Protein antibody (Biotin) (orb453363)

Description	Maltose Binding Protein antibody (Biotin)
Species/Host	Mouse
Reactivity	E. coli, Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant MBP Protein of Escherichia coli (strain K12)
Target	MBP tag
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	POAEX9
Dilution Range	WB=1:1000-5000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.