

Product Datasheet

Scavenger Receptor BI + BII antibody (Biotin) (orb455286)

Description	Scavenger Receptor BI + BII antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Scavenger Receptor BI + BII
Target	Scavenger Receptor BI + BII
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WTV0
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.