

Product Datasheet

PID1 antibody (FITC) (orb4554)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal to PID1. Increases proliferation of preadipocytes without affecting adipocytic dif... |
| Species/Host | Rabbit |
| Conjugation | FITC |
| Tested Applications | IF |
| Immunogen | KLH conjugated synthetic peptide derived from human PID1 |
| Target | PID1 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Lyophilized or Liquid |
| Concentration | - 100 µg (in 100 µl): 1 mg/ml |
| Storage | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| MW | 28kDa |
| Uniprot ID | Q7Z2X4 |
| Expiration Date | 12 months from date of receipt. |